


ApneaLink Air Patient Instructions

Video Instructions: <https://www.youtube.com/watch?v=MNRxmhdLjW0&feature=youtu.be>

Equipment Pieces:

- ApneaLink Air Device (#1)
- Effort sensor (#2)
- Belt (#3)
- Oximeter belt clip (optional - #4)
- Oximeter adapter cable (#5 & #6)
- Finger Pulse Oximeter (#7)
- Nasal Cannula (#8) (new)
- ResMed Protective Bag (not shown)
- Rechargeable Batteries (🔋 **MUST RETURN WITH EQUIPMENT - A CHARGE OF \$10.00 WILL APPLY TO MISSING BATTERIES**)



Change device batteries after each night.

